



Rheomat R 190

The viscometer for
product control in
laboratory and field



proRheo
rheological answers!

COMPATIBLE

New USB port for software use and measurements directly from PC.

DURABLE

The housing of the R190 is made of solid aluminium and offers the usual proRheo quality. It will meet all the challenges of everyday work.

FLEXIBLE MEASUREMENT MODES

• Single point measurement:

Enter any shear rate and read the measurement result immediately. If your R190 is connected to a PC the measurement results are recorded continuously. Changes in shear stress and viscosity are also recorded as a function of time. (Excellent for pot life determination).

• Step Programs:

In the automatic mode you can choose between 10 measurement programs, recording 8 points at different shear rates. The rotational speed is at first increased from low to high and then back down again. The resulting points give a flow curve which can be recorded on your PC or printer. Two of the step programs are predefined and a further 8 are available for customised testing, these can be saved and retrieved by their step program number.



NEW DESIGN

The R190 features a modern design with a bright LED color touch display with swivel function.

OVERVIEW

The following values are displayed and continuously updated:

- Temperature
- Measurement System Number
- Torque and shear rate
- Rotational speed and shear stress
- Viscosity

SAFE TO HANDLE

A clear dialog guides you through the necessary input options. The buttons marked „Manual“, „Automatic“, „Printer“ and „PC“ start the respective functions directly.

FIELD-PROVEN

The R190 uses an integrated grip built right into the housing for ease of use.

MEMORY

The integrated measured value memory is powered by a separate lithium battery and saves your measured values.

PRACTICAL

The R190 is supplied with all accessories in a case and is thus quickly ready for use anywhere.



R190 DIMENSIONS

Weight: 2,55 kg

Dimension: 150 x 345 x 130 (W x H x D/mm)

INSTRUMENT OPERATIONAL INFORMATION

The equipment may be stored and operated in an environment from -20 to 60 °C.

VOLTAGE

with power supply:

100 to 250 V AC with 50/60 Hz, 1.0 A 

INTERFACES

USB for direct PC connection.

TORQUE

0,1 to 25 mNm +/- 0,01 mNm

ROTATIONAL SPEED

5 to 1000 rpm +/- 1 rpm

MEASURING SYSTEMS

11 predefined measuring systems

99 programmable measuring systems

MEASUREMENT RANGE

Viscosity: 0,002 Pas to 10.000 Pas

according to measurement systems.

Share range: 0,8 s⁻¹ to 3.000 s⁻¹

MEASUREMENT OF TEMPERATURE BY PT 100

TEMPERATURE OF SAMPLE

-9,9 to 99,9 °C +/- 0,1 °C

100 to 120 °C: +/- 1,0 °C

MEASUREMENT PROGRAMS

- 8 measurement points at different shear rates.
- 2 predefined test setups.
- 8 programmable test setups, minimum and maximum shear rates, statistical analysis available.

SOFTWARE

Different operational and analytical programs are available. Special print-out programs for ASCII or direct access to Excel.

System requirements for software RHESY:

- IBM PC or compatible PC from Pentium
- 166 MHz, 64 MB main memory (RAM)
- CD-ROM drive
- 1 free serial or USB interface
- Windows operating system

	Measurement systems	Measurement tube Ø mm	Measurement bob Ø mm	Viscosity (Pas)min.	Viscosity (Pas)max.	filling volume (ml)
DIN 53018/ DIN 53019	11 22 33	32,54 26,03 15,18	30 24 14	0,005 0,010 0,050	19 38 191	ca. 24 ca. 16 ca. 9
Relative systems	19 12 13 23 14	32,54 32,54 32,54 26,03 32,54	31,5 24 14 14 14	0,002 0,027 0,210 0,240 0,545	7 104 800 906 2.080	ca. 20 ca. 18 ca. 26 ca. 18 ca. 26
Special relative systems	71 71 73 74 75			0,003 0,027 0,160 0,665 2,580	10 104 605 2.530 9.800	
ISO 2555	61 62 63 64 65 66 67			0,007 0,028 0,070 0,139 0,278 0,696 2,783	26 106 264 529 1.057 2.643 10.574	